ENZYWELL

EPSTEIN BARR EA IgG

REF 91058 (96 tests)

Prodotto da/Manufactured by/Fabricado por:
DIESSE Diagnostica Senese
Via delle Rose 10
53035 Monteriggioni (Siena) - Italy

INDICE / INDEX / INDICE / CONTEÚDO

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17. BIBLIOGRAFIA

- Gorgievski H.M., Hinderer W. et al. Serodiagnosis of infectious mononucleosis by using recombinant Epstein-Barr virus antigens and enzyme-linked immunosorbent assay technology. *J. Clin. Microbiol.* 28: 2305 (1990).
- Middeldorp, J.M. and Herbrink, P. Epstein-Barr virus specific marker molecules for early diagnosis of infectious mononucleosis. *J. Virol. Methods* 21: 133-146 (1988).
- Luka, J., Chase, R.C. and Pearson G.R. A sensitive enzyme linked immunosorbent assay (ELISA) against the major EBV-associated antigens. I. Correlation between ELISA and immunofluorescence titers using purified antigens. *J. Immunol. Methods* 67: 145.156 (1984).
- Motz, M., Fan J., Seibl R., and Wolf H. Expression of the Epstein Barr Virus 138 KDa early protein in *E. coli* for the use as antigen in diagnostic tests. *Gene* 42: 303-312 (1986).

DIESSE

INSTRUCTIONS FOR USE

ENZYWELL
EPSTEIN BARR
EA IgG

REF 91058

(English)

1. INTENDED USE

IMMUNOENZYMATIC METHOD FOR THE QUALITATIVE DETERMINATION OF IgG-CLASS ANTIBODIES TO EPSTEIN BARR VIRUS EARLY ANTIGEN IN HUMAN SERUM, AS AN AID IN THE DIAGNOSIS OF EPSTEIN BARR VIRUS.

2. SUMMARY AND EXPLANATION OF THE TEST

The Epstein Barr Early Antigens (EBV-EA) are viral proteins which are produced in the early stages of lythic infection. In the diagnosis of infectious mononucleosis (IM), indirect immunofluorescence has been adopted for a long time for the determination of serum antibodies against EBV-EA, in association with the detection of IgG and IgM against the capsid antigens (VCA), and IgG against a nuclear antigen (EBNA-1). Recently it has been shown that the immunoenzymatic determination of serum IgG and IgM against the two early proteins, p138 and p54, shows good correlation with the indirect immunofluorescence method and can be utilised, not only as a complementary method in the determination of anti-VCA antibodies, but also as a substitute for these methods in the diagnosis of IM.

The DIESSE Kit for the immunoenzymatic determination of IgG and IgM anti-EBV-EA uses recombinant polypeptides produced in E. coli containing epitopes of the proteins p138 and p54.

3. PRINCIPLE OF THE TEST

The test is based on the ELISA technique (Enzyme linked Immunosorbent Assay).

The recombinant antigens are bound to the solid phase. Following incubation with dilute human serum, the specific immunoglobulins are bound to the antigen.

After washings to eliminate the proteins which have not reacted, incubation with the conjugate, composed of human IgG monoclonal antibodies labelled with peroxidase, is performed.

The unbound conjugate is eliminated and the peroxidase substrate is added.

The colour which develops is proportional to the concentration of specific antibodies present in the serum sample.

4. KIT CONTENTS AND REAGENT PREPARATION

- Reagents are sufficient for 96 determinations.

- **Bring reagents to room temperature before use.**

MT PLATE MICROPLATE 1 x 96 wells coated with recombinant EA proteins.

Use: open the package at the opposite end from the code (I, followed by the lot number) which is useful for identification purposes, remove the support and wells to be used from the foil package, and place the unused wells in the polythene bag with the silica gel, expel the air and seal by pressing the closure.

CONJ CONJUGATE 1 x16 mL

Contents: anti human IgG monoclonal antibodies labelled with Peroxidase, in phosphate buffer solution containing phenol 0.05% and Bronidox 0.02%.

Preparation: ready for use.

CONTROL + POSITIVE CONTROL 1 x 1.6 mL

Contents: Diluted human serum at known concentration of anti-EA IgG, in Phosphate buffer 0.01 mol/L with BSA 1% and sodium azide 0.09%, liquid, ready for use without further dilution.

Colour: the colour is proportional to the relative antibody titer.

CONTROL CUT-OFF CUT OFF CONTROL 1 x 2 mL

Contents: Diluted human serum at known concentration of anti-EA IgG, in Phosphate buffer 0.01 mol/L with BSA 1% and sodium azide 0.09%, liquid, ready for use without further dilution.

Colour: the colour is proportional to the relative antibody titer.

CONTROL IgG - IgG NEGATIVE CONTROL (PF93910) 1 x 1.6 mL **INTERCHANGEABLE BETWEEN LOTS OF THE EBV PANEL, REF 91055 VCA IgG, 91056 VCA IgM, 91057 EBNA IgG, 91058 EA IgG, 91059 EA IgM**

Contents: Diluted human serum in Phosphate buffer 0.01 mol/L with BSA 1% and sodium azide 0.09%, liquid, ready for use without further dilution.

WASH BUF 10x WASH BUFFER 10X (PF93603) 1 x 100 mL **INTERCHANGEABLE BETWEEN LOTS**

Contents: Phosphate buffered saline, concentrated 10 times; contains Brij 0.5%.

Preparation: dilute the required volume 1:10 with distilled water in order to obtain the washing buffer ready for use. If crystals are present, they should be dissolved at 37°C before dilution.

SAMP DIL 10 DILUENT 10 (PF93621) 1 x 100 mL **INTERCHANGEABLE BETWEEN LOTS**

To be used to dilute samples.

Contents: PBS containing proteins 10% w/v and sodium azide 0,09% containing methyl orange as dye.

SORBENT G EBV DILUENT (PF91083) 1 x 7 mL **INTERCHANGEABLE BETWEEN LOTS**

To be placed in the wells to neutralize non specific reactions.

Contents: Phosphate buffered saline (PBS) containing sodium azide 0.09% and bovine serum 10% w/v. Ready for use.

SUBS TMB SUBSTRATE (PF93619) 1x 12 mL Ready for use **INTERCHANGEABLE BETWEEN LOTS**

Contents: Tetramethylbenzidine 0.26 mg/mL and hydrogen peroxide 0.01% stabilised in citrate buffer 0.05 mol/L (pH 3.8).

H₂SO₄ 0.3 M STOP SOLUTION (PF93602) 1x16 mL **INTERCHANGEABLE BETWEEN LOTS**

H₂SO₄ 0.3 mol/L, in solution ready for use.

ADHESIVE FILMS (2).

POLYTHENE BAG (1).

MATERIALS REQUIRED BUT NOT PROVIDED.

- Incubator at 37°C
- Microplate reader (wave length 450 nm or 450/620 nm, with linearity up to OD >= 2000)
- Microplate washer (preferable) able to dispense volumes in the range 225-375 µL
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 10, 100, 1000 µL solution
- Disposable gloves
- Timer
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials
- Absorbent tissue.

5. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C.

The expiry date is printed on each component and on the box label.

Reagents have a limited stability after opening and/or preparation

REAGENT	CONDITIONS
Microplate	5 weeks at 2/8°C, polythene bag
Controls	5 weeks at 2/8°C
Conjugate	5 weeks at 2/8°C
Substrate	up to the expiry date at 2/8°C, 1 week at 15-30°C; store in the dark
Sample Diluent	up to the expiry date at 2/8°C
Wash Buffer	2 weeks at 2/8°C, 5 days at 15/30°C.
Sorbent G	5 weeks at 2/8°C
Stop Solution	up to the expiry date at 2/8°C

6. PRECAUTIONS **FOR IN VITRO DIAGNOSTIC USE ONLY.**

Caution:

This kit contains materials of human origin which have been tested and gave a negative response by FDA-approved methods for the presence of HbsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: serum samples and reagents once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth. Wear disposable gloves and eye protection while handling specimens and performing the assay. Wash hands thoroughly when finished.
2. The following reagents contain low concentrations of harmful or irritant substances:
 - a) The Wash Buffer contains detergents
 - b) The conjugate contains phenol
 - c) The substrate is acid
 - d) The controls contain 0.09% Sodium Azide which can react with lead and copper in plumbing forming highly explosive deposits of metal azides.

If any of the reagents come into contact with the skin or eyes, wash the area extensively with water.
3. Non-disposable apparatus should be sterilized after use. The preferred method is to autoclave for 1 h at 121°C; disposables should be autoclaved or incinerated.
4. Sulphuric acid required for the Stop Solution and hydrochloric acid used for washing glassware are corrosive and should be handled with appropriate care. If they come into contact with the skin or eyes, wash thoroughly with water.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1.0%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1.0% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical precautions

1. Allow all reagents and samples to come to room temperature (18-30°C) before use. Immediately after use return reagents to the recommended storage temperature. **It is important to work at the correct temperature. Check that the thermostat does not go below 35°C or over 39°C.** Open the envelope containing the strips after at least ½ hr at room temperature.
2. Do not use the reagents beyond the stated expiry date. Microbiological contamination of reagents must be avoided as this may reduce the life of the product and cause erroneous results.
3. Do not modify the Test Procedure or substitute reagents from other manufacturers or other lots unless the reagent is stipulated as interchangeable. Do not reduce any of the recommended incubation times.
4. Any glassware to be used with the reagents should be thoroughly washed with 2M hydrochloric acid and then rinsed with distilled water or high quality deionized water.
5. Do not expose reagents to strong light or hypochlorite fumes during storage or during incubation steps.
6. Do not allow wells to become dry during the assay procedure.
7. Care must be taken not to cross-contaminate reagents. It is important that pipettes are dedicated for exclusive use with the various reagents.
8. Care should be taken to avoid touching or splashing the rim of the well with conjugate. Do not "blow-out" from microplates.
9. Enzyme immunoassays can occasionally exhibit an "edge effect" which must be minimized by increasing the humidity during incubation steps. Plates must be covered with their covers and incubated at 37°C either in a water bath with a rack or float to support the plates if necessary, or in an incubator. Alternatively, plates can be incubated in an approved analyzer. See the appropriate operating manual for further details. CO₂ incubators must not be used.
10. Ensure that the bottom of the plate is clean and dry, and that no bubbles are present on the surface of the liquid before reading the plate.
11. Use of highly hemolyzed samples, incompletely clotted sera, or samples with microbial contamination may give rise to erroneous results.

12. Use of the kit on automatic instruments must be validated by the user.
13. For each instrument used, read the manufacturer's instructions manual carefully to obtain additional information on the following points:
- installation and particular requisites
 - operating principles, instructions, precautions and risks
 - manufacturer's specifications and instrument performance
 - servicing and maintenance.

7. TYPE AND STORAGE OF SAMPLE

The sample is composed of serum obtained from blood collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice. The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at -20°C, and can be thawed a maximum of 3 times. Do not keep the samples in auto-defrosting freezers. Defrosted samples must be carefully mixed before performing the test. Heat inactivation can lead to erroneous results. The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

Strongly lipemic, icteric or contaminated samples should be avoided. **The test is not applicable to human plasma.**

8. TEST PROCEDURE

Manual Technique

- Prepare the required number of strips.
- Prepare the washing buffer by diluting the Wash Buffer 10x (100 mL + 900 mL H₂O).

Dilute samples 1:26 distributing 40 µL of serum into 1 mL of diluent. Leave one well for the blank, performed using 100 µL of the substrate mixture. Dispense 50 µL of Sorbent G and 50 µL of each diluted sample per well (duplicate testing is recommended). Place UNDILUTED controls in a strip (100 µL in each well without the Sorbent G). The minimum requisite is 1 negative control, 2 cut-off and 1 positive control.

Wells are covered with protective film and incubated for 45 minutes at 37°C. After washing four times for 30 seconds (300 µL), add 100 µL of the conjugate to each well and incubate again for 45 minutes at 37°C, covering the wells with the protective film. The plate is washed again 4 times, as described above. Finally, the substrate is distributed, 100 µL/well.

After 15 minutes at room temperature the enzymatic reaction is stopped with 100 µL of Stop Solution.

The adsorbance (O.D.) is read at 450 nm or at 450/620 nm within 30 min.

9. SCHEME OF TEST PROCEDURE

- | | |
|--------|---|
| STEP 1 | Place 50 µL of Sorbent G, 50 µL of diluted sample and 100 µL controls in the wells of the strips. |
| | - |
| | Incubate for 45 min. at 37°C |
| | - |
| | Wash 4 times (300 µL) |
| | - |
| STEP 2 | Add 100 µL of conjugate to each well |
| | - |
| | Incubate for 45 min. at 37°C |
| | - |
| | Wash 4 times (300 µL) |
| | - |
| STEP 3 | Add 100 µL of Substrate to each well |
| | - |
| | Incubate for 15 min. at R.T. |
| | - |
| STEP 4 | Add 100 µL of Stop Solution |
| | - |
| | Read absorbance at 450 nm within 30 min |

10. TEST VALIDATION

Subtract the value of the blank (≤ 0.150) from all the other readings. The Positive control must have an O.D. at least 1.5 times that of the Cut-Off serum. The ratio between Negative Control and Cut-off must be ≤ 0.6 . The O.D. of the Cut-off must be ≥ 0.2 at 450 nm and ≥ 0.16 at 450/620 nm.

11. INTERPRETATION OF THE RESULTS

Qualitative results

If the adsorbance of the sample is higher than that of the Cut-off, the sample is positive for the presence of specific IgG. Calculate the ratio between the O.D. value of the sample and that of the Cut-off (INDEX). The sample is considered:

Positive: if the ratio is > 1.2 .

Doubtful: $\pm 20\%$ of the Cut-off.

Negative: if the ratio is < 0.8 .

If the result is doubtful, repeat the test. If it remains doubtful, collect a new serum sample.

Quantitative results

Arbitrary Units (AU)

Positive results can be expressed in AU by applying the following formula:

$$\text{AU} = \frac{\text{O.D. sample}}{\text{O.D. Cut-Off}}$$

This ratio is generally considered as an "index" of positivity.

12. LIMITATIONS OF THE PROCEDURE

The results obtained can only be used as a guide in the diagnosis and must always be evaluated together with results of other diagnostic procedures.

There is the possibility of cross-reactions with samples containing anti-E. coli antibodies.

Performance characteristics have not been determined for a visual interpretation, nor have they been evaluated in patients affected by nasopharyngeal carcinoma, Burkitt's lymphoma or other lymphadenopathies associated with EBV, apart from EBV-related mononucleosis.

As anti-EA IgG antibodies are present in normal convalescent sera, the diagnosis cannot be established on the basis of a single parameter. An accurate interpretation of an EBV infection must be based on the results of VCA IgG, VCA IgM, EBNA IgG, EA IgG ed EA IgM results, in addition to heterophyle antibodies.

The test must be performed on serum. The use of blood or plasma has not been established.

13. ANALYTICAL SPECIFICITY

10 samples which were negative for anti-EA IgG and which contained IgG anti-CMV, anti-Herpes Simplex and anti-Varicella, were tested. In no case did the presence of these antibodies interfere with the Enzywell EA IgG kit.

14. DIAGNOSTIC SENSITIVITY AND SPECIFICITY

118 samples were analysed, 58 of which were negative and 60 positive. The comparative method (another commercial ELISA kit) gave 67 and 51 samples, respectively; those which were not in agreement were tested with the other serological parameters typical of Epstein Barr Virus infection, and 5 showed the typical pattern of an acute EBV infection in progress, confirming the positive results.

The results are the following:

		REFERENCE	
		+	-
DIESSE	+	56	4
	-	1	57

The ENZYWELL method has a sensitivity of 98.2% and specificity of 93.4%.

15. PRECISION

“Within run” precision between 3 different lots

<i>Cut-off n=15</i>	<i>Lot 014</i>	<i>Lot 015</i>	<i>Lot 016</i>
<i>O.D.</i>	0.452	0.569	0.522

CV%	8	6	10
-----	---	---	----

“Between run” and “between lots” precision

Samples	INDEX			Average	CV%
	Lot 014	Lot 015	Lot 016		
1	0.3	0.2	0.2	0.2	25
2	1.9	1.7	1.6	1.7	9
3	2.8	2.8	2.6	2.7	4

16. TROUBLE SHOOTING GUIDE

PROBLEM	POSSIBLE SOURCE	TEST OR ACTION
Invalid run (all negative)	One or more reagents not added or added in wrong sequence	Recheck procedure Check for unused solutions. Repeat test.
	Unreactive plate	Check the code on the package containing the plate (see package insert point 4 for correct code).
		Check for moisture in unused plate. (Silica gel dessiccant must be pale yellow). Repeat test
Invalid run (all positive)	Contamination of substrate	Take new aliquot of substrate.
	Inadequate washing	Ensure that wash apparatus works well
Poor precision	Incomplete washing of wells	Ensure that wash apparatus works well
	Inadequate aspiration of wells	Ensure that wash apparatus works well
	Pipetting error	Check pipette function
	Reagent addition too slow	Avoid drying of the plate after washing step. Add reagents immediately
	Presence of bubbles	Avoid air bubbles during pipetting.
	Optical pathway not clean	Check instrument light source and detector for dirt. Wipe bottom of plate with soft tissue.
Inadequate Color development	Incorrect incubation times or temperature	Check for temperature control and time monitoring
		Adhere to recommended instruction for use.
	Inadequate volume of substrate added to the plate	Check pipette function.

17. REFERENCES

- Gorgievski H.M., Hinderer W. et al. Serodiagnosis of infectious mononucleosis by using recombinant Epstein-Barr virus antigens and enzyme-linked immunosorbent assay technology. J. Clin. Microbiol. 28: 2305 (1990).
- Middeldorp, J.M. and Herbrink, P. Epstein-Barr virus specific marker molecules for early diagnosis of infectious mononucleosis. J. Virol. Methods 21: 133-146 (1988).
- Luka, J., Chase, R.C. and Pearson G.R. A sensitive enzyme linked immunosorbent assay (ELISA) against the major EBV-associated antigens. I. Correlation between ELISA and immunofluorescence titers using purified antigens. J. Immunol. Methods 67: 145.156 (1984).
- Motz, M., Fan J., Seibl R., and Wolf H. Expression of the Epstein Barr Virus 138 KDa early protein in E. coli for the use as antigen in diagnostic tests. Gene 42: 303-312 (1986).

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1. Gorgievski H.M., Hinderer W. et al. Serodiagnosis of infectious mononucleosis by using recombinant Epstein-Barr virus antigens and enzyme-linked immunosorbent assay technology. *J. Clin. Microbiol.* 28: 2305 (1990).
2. Middeldorp, J.M. and Herbrink, P. Epstein-Barr virus specific marker molecules for early diagnosis of infectious mononucleosis. *J. Virol. Methods* 21: 133-146 (1988).
3. Luka, J., Chase, R.C. and Pearson G.R. A sensitive enzyme linked immunosorbent assay (ELISA) against the major EBV-associated antigens. I. Correlation between ELISA and immunofluorescence titers using purified antigens. *J. Immunol. Methods* 67: 145.156 (1984).
4. Motz, M., Fan J., Seibl R., and Wolf H. Expression of the Epstein Barr Virus 138 KDa early protein in *E. coli* for the use as antigen in diagnostic tests. *Gene* 42: 303-312 (1986).

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

17. REFERÊNCIAS

1. Gorgievski H.M., Hinderer W. et al. Serodiagnosis of infectious mononucleosis by using recombinant Epstein-Barr virus antigens and enzyme-linked immunosorbent assay technology. *J. Clin. Microbiol.* 28: 2305 (1990).
2. Middeldorp, J.M. and Herbrink, P. Epstein-Barr virus specific marker molecules for early diagnosis of infectious mononucleosis. *J. Virol. Methods* 21: 133-146 (1988).
3. Luka, J., Chase, R.C. and Pearson G.R. A sensitive enzyme linked immunosorbent assay (ELISA) against the major EBV-associated antigens. I. Correlation between ELISA and immunofluorescence titers using purified antigens. *J. Immunol. Methods* 67: 145-156 (1984).
4. Motz, M., Fan J., Seibl R., and Wolf H. Expression of the Epstein Barr Virus 138 KDa early protein in *E. coli* for the use as antigen in diagnostic tests. *Gene* 42: 303-312 (1986).

	EN ES IT	Date of manufacture Fecha de fabricación Data di fabbricazione	FR GR PT	Date de fabrication Ημερομηνία Παραγωγής Data de fabrico
	EN ES IT	Use By Fecha de caducidad Utilizzare entro	FR GR PT	Utiliser jusque Ημερομηνία λήξης Prazo de validade
	EN ES IT	Caution, consult accompanying documents Atención, ver instrucciones de uso Attenzione, vedere le istruzioni per l'uso	FR GR PT	Attention voir notice d'instructions Προειδοποίηση, συμβουλευτείτε τα συνοδά έντυπα Atenção, consulte a documentação incluída
	EN ES IT	Manufacturer Fabricante Fabbricante	FR GR PT	Fabricant Κατασκευαστής Fabricante
	EN ES IT	Contains sufficient for <n> tests Contenido suficiente para <n> ensayos Contenuto sufficiente per "n" saggi	FR GR PT	Contenu suffisant pour "n" tests Περιεχόμενο επαρκές για «n» εξετάσεις Conteúdo suficiente para "n" ensaios
	EN ES IT	Temperature limitation Límite de temperatura Limiti di temperatura	FR GR PT	Limites de température Περιορισμοί θερμοκρασίας Limites de temperatura
	EN ES IT	Consult Instructions for Use Consulte las instrucciones de uso Consultare le istruzioni per l'uso	FR GR PT	Consulter les instructions d'utilisation Συμβουλευτείτε τις οδηγίες χρήσης Consulte as instruções de utilização
	EN ES IT	Biological risks Riesgo biológico Rischio biologico	FR GR PT	Risques biologiques Βιολογικοί κίνδυνοι Risco biológico
	EN ES IT	Catalogue number Número de catálogo Numero di catalogo	FR GR PT	Référence du catalogue Αριθμός καταλόγου Referência de catálogo
	EN ES IT	In Vitro Diagnostic Medical Device Producto sanitario para diagnóstico in vitro Dispositivo medico-diagnostico in vitro	FR GR PT	Dispositif médical de diagnostic in vitro In Vitro Διαγνωστικό Ιατροτεχνολογικό προϊόν Dispositivo médico para diagnóstico in vitro
	EN ES IT	Batch code Código de lote Codice del lotto	FR GR PT	Code du lot Αριθμός Παρτίδας Código do lote

DIESSE Diagnostica Senese
Via delle Rose 10
53035 Monteriggioni (Siena) Italy
Tel. 0577-587111