*ENZYWELL*

**CYTOMEGALOVIRUS IgM**

**REF** 91011 (96 tests)

Prodotto da/Manufactured by/Fabricado por:  
DIESSE Diagnostica Senese  
Via delle Rose 10  
53035 Monteriggioni (Siena) - Italy

**INDICE / INDEX / INDICE / CONTEÚDO**

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## **17. BIBLIOGRAFIA**

1. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. H.O. Kangro: Cytomegalovirus serology: does it give the answer? Serodiagnosis and Immunotherapy 1: 91 (1987).
3. Grint P.C.A. et al.: Screening tests for antibodies to cytomegalovirus: an evaluation of five commercial products. J. Clin. Pathol. 38: 1059 (1985).
4. Van Loon A.M. et al.: Direct enzyme-linked immunosorbent assay that uses peroxidase-labelled antigen for determination of immunoglobulin M antibody to cytomegalovirus. J. Clin. Microbiol. 13: 416
5. M. Musiani et al.: Rapid detection of antibodies against cytomegalovirus induced immediate early and early antigens by an enzyme linked immunosorbent assay. J. Clin. Pathol. 37: 122 (1984).
6. F. de Ory et al.: Serological diagnosis of cytomegalovirus infections: comparison of six commercial methods of ELISA. Serodiagnosis and Immunotherapy 2: 423 (1988).
7. P. Dal Monte et al.: Cytomegalovirus umano. Diagnosis 2: 67 (1990).
8. R. Ziegelmaier et al.: ELISA. la Ricerca Clin. Lab. 10: 83 (1980).

**RIEDEL  
PHARMACIA**  
**INSTRUCTIONS FOR USE**

**ENZYWELL  
CYTOMEGALOVIRUS IgM**

REF 91011

(English)

**1. INTENDED USE**

**IMMUNOENZYMATIC CAPTURE METHOD FOR THE QUALITATIVE DETERMINATION OF IgM-CLASS ANTIBODIES TO CYTOMEGALOVIRUS IN HUMAN SERUM, AS AN AID IN THE DIAGNOSIS OF CYTOMEGALOVIRUS INFECTION.**

**2. SUMMARY AND EXPLANATION OF THE TEST**

Cytomegalovirus is a herpes virus transmitted by close human contact. No symptoms of infection are apparent in the majority of cases. However, the virus is very dangerous and may be fatal in immunodepressed patients. Serum-negative female patients who become infected during pregnancy may transmit the disease to the fetus. In 95% of cases this occurs without symptoms, but some neonates may present jaundice, hepato-splenomegaly and retarded psycho-motorial development. For this reason it is of great importance to determine the immunitary state of the patient before the start of pregnancy, if possible, and check for serum conversion.

The assay of specific IgM is of great importance in the diagnosis of primary infection.

**3. PRINCIPLE OF THE TEST**

The test for the assay of Cytomegalovirus IgM is based on the principle of the capture of these immunoglobulins by anti-human IgM monoclonal antibodies found on the solid phase. A subsequent incubation with cytomegalovirus antigen in a complex with monoclonal antibodies conjugated to horse radish peroxidase selects the IgM antibodies specific for the antigen and is revealed by the addition of the peroxidase substrate. When the enzymatic reaction is stopped by the addition of a sulphuric acid solution, a yellow colouring forms. The colour, which is proportional to the amount of specific antibodies present in the sample, can be read in an ELISA microplate reader.

**4. KIT CONTENTS AND REAGENT PREPARATION**

- Reagents are sufficient for 96 determinations.

**Bring to room temperature before use.**

**MT PLATE** MICROPLATE 12 x 8 wells coated with anti-human IgM monoclonal antibodies.

Use: open the package on the opposite side from the code (M followed by the lot number) which is useful for identification purposes, remove the support and strips to be used from the foil package, and place the unused strips in the polythene bag with the silica gel, expel the air and seal by pressing the closure.

**CONTROL +** POSITIVE CONTROL 1 x 1.6 mL

Contents: Diluted human serum in Phosphate buffer 0.01 mol/L with BSA 1% and sodium azide 0.09%, liquid, ready for use without further dilution.

Colour: the colour of the controls is proportional to the relative antibody titer.

**CONTROL CUT-OFF** CUT OFF CONTROL 1 x 2.5 mL

Contents: Diluted human serum in Phosphate buffer 0.01 mol/L with BSA 1% and sodium azide 0.09%, liquid, ready for use without further dilution.

Colour: the colour of the controls is proportional to the relative antibody titer.

**Ag** ANTIGEN Freeze-dried powder x 6 vials

Contents: Partially purified Cytomegalovirus, inactivated by treatment with Beta-propiolactone, in Phosphate buffer containing mouse ascitic fluid and lactose.

Preparation: reconstitute with the conjugate volume shown on the label, mixing by inversion.

**CONJ** CONJUGATE 1 x 18 mL

**Contents:** monoclonal antibodies labelled with peroxidase, in phosphate buffer with phenol 0.05% and Bronidox 0.02%.

**Preparation:** ready for use.

The immunocomplex should be prepared about 45 min. before use.

**CONTROL IgM** – IgM NEGATIVE CONTROL (PF93900) 1 x 1.6 mL **INTERCHANGEABLE BETWEEN LOTS**

**Contents:** Diluted human serum in Phosphate buffer 0.01 mol/L with BSA 1% and sodium azide 0.09%, liquid, ready for use without further dilution.

**WASH BUF 10x** WASH BUFFER 10X (PF93603) 1 x 100 mL **INTERCHANGEABLE BETWEEN LOTS**

**Contents:** Phosphate buffered saline, concentrated 10 times; contains Brij 0.5% .

**Preparation:** dilute the required volume 1:10 with distilled water in order to obtain the washing buffer ready for use. If crystals are present, they should be dissolved at 37°C before dilution.

**SAMP DIL 2** DILUENT 2 (PF93611) 1 x 100 mL **INTERCHANGEABLE BETWEEN LOTS**

To be used to dilute samples.

**Contents:** Proteic solution in phosphate buffer with sodium azide 0,09% containing methyl orange as dye.

**SUBS TMB** SUBSTRATE (PF93619). 1x12 mL. Ready for use. **INTERCHANGEABLE BETWEEN LOTS**

**Contents:** Tetramethylbenzidine 0.26 mg/mL and hydrogen peroxide 0.01% stabilised in citrate buffer 0.05 mol/L (pH 3.8).

**H<sub>2</sub>SO<sub>4</sub> 0.3 M** STOP SOLUTION (PF93602) 1x16 mL **INTERCHANGEABLE BETWEEN LOTS**

H<sub>2</sub>SO<sub>4</sub> 0.3 mol/L, in solution ready for use.

ADHESIVE FILMS (2).

POLYTHENE BAG (1).

**MATERIALS REQUIRED BUT NOT PROVIDED.**

- Incubator at 37°C
- Microplate reader (wave length 450 or 450/620 nm, with linearity up to OD >= 2000)
- Microplate washer (preferable) able to dispense volumes in the range 225-375 µL
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 10, 100, 1000 µl solution
- Disposable gloves
- Timer
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials
- Absorbent tissue.

**5. STORAGE AND STABILITY OF REAGENTS**

Reagents must be stored at 2/8°C.

The expiry date is printed on each component and on the box label.

**Reagents have a limited stability after opening and/or preparation**

REAGENT	CONDITIONS
Microplate	5 weeks at 2/8°C, polythene bag
Controls	5 weeks at 2/8°C
Conjugate	5 weeks at 2/8°C
Reconstituted antigen	15 days at 2/8°C,
Substrate	until the expiry date at 2/8°C, 1 week at 15-30°C in the dark
Sample Diluent	up to the expiry date at 2/8°C
Wash Buffer	2 weeks at 2/8°C, 5 days at 15/30°C.
Stop Solution	until the expiry date at 2/8°C

**6. PRECAUTIONS**

FOR IN VITRO DIAGNOSTIC USE ONLY.

**Caution:** This kit contains materials of human origin which have been tested and gave a negative response by FDA-approved methods for the presence of HbsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

**Waste disposal:** serum samples and reagents once used must be treated as infectious residuals and eliminated according to law.

#### Health and Safety Information

1. Do not pipette by mouth. Wear disposable gloves and eye protection while handling specimens and performing the assay. Wash hands thoroughly when finished.
2. The following reagents contain low concentrations of harmful or irritant substances:
  - a) The Wash Buffer contains detergents
  - b) The conjugate contains phenol
  - c) The substrate is acid
  - d) The controls contain 0.09% Sodium Azide which can react with lead and copper in plumbing forming highly explosive deposits of metal azides.

If any of the reagents come into contact with the skin or eyes, wash the area extensively with water.
3. Non-disposable apparatus should be sterilized after use. The preferred method is to autoclave for 1 h at 121°C; disposables should be autoclaved or incinerated.
4. Sulphuric acid required for the Stop Solution and hydrochloric acid used for washing glassware are corrosive and should be handled with appropriate care. If they come into contact with the skin or eyes, wash thoroughly with water.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1.0%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1.0% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

#### Analytical precautions

1. Allow all reagents and samples to come to room temperature (18-30°C) before use. Immediately after use return reagents to the recommended storage temperature. **It is important to work at the correct temperature. Check that the thermostat does not go below 35°C or over 39°C.** Open the envelope containing the strips after at least ½ hr at room temperature.
2. Do not use the reagents beyond the stated expiry date. Microbiological contamination of reagents must be avoided as this may reduce the life of the product and cause erroneous results.
3. Do not modify the Test Procedure or substitute reagents from other manufacturers or other lots unless the reagent is stipulated as interchangeable. Do not reduce any of the recommended incubation times.
4. Any glassware to be used with the reagents should be thoroughly washed with 2M hydrochloric acid and then rinsed with distilled water or high quality deionized water.
5. Avoid the use of self-defrosting freezers for the storage of samples.
6. Do not expose reagents to strong light or hypochlorite fumes during storage or during incubation steps.
7. Do not allow wells to become dry during the assay procedure.
8. Care must be taken not to cross-contaminate reagents. It is important that pipettes are dedicated for exclusive use with the various reagents.
9. Care should be taken to avoid touching or splashing the rim of the well with conjugate. Do not "blow-out" from microplates.
10. Enzyme immunoassays can occasionally exhibit an "edge effect" which must be minimized by increasing the humidity during incubation steps. Plates must be covered with their covers and incubated at 37°C either in a water bath with a rack or float to support the plates if necessary, or in an incubator. Alternatively, plates can be incubated in an approved analyzer. See the appropriate operating manual for further details. CO<sub>2</sub> incubators must not be used.
11. Ensure that the bottom of the plate is clean and dry, and that no bubbles are present on the surface of the liquid before reading the plate.
12. Use of highly hemolyzed samples, incompletely clotted sera, or samples with microbial contamination may give rise to erroneous results.
13. The use of the kit on automatic instruments must be validated by the user.
14. For each instrument used, read the manufacturer's instructions manual carefully to obtain additional information on the following points:

- installation and particular requisites
- operating principles, instructions, precautions and risks
- manufacturer's specifications and instrument performance
- servicing and maintenance.

## **7. TYPE AND STORAGE OF SAMPLE**

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice. The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at -20°C, and can be thawed a maximum of 3 times. Do not keep the samples in auto-defrosting freezers. Defrosted samples must be carefully mixed before performing the test. Heat inactivation can lead to erroneous results. The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

Strongly lipemic, icteric or contaminated samples should be avoided.

**The test is not applicable to human plasma.**

## **8. TEST PROCEDURE**

### **Manual Technique**

- Prepare the required number of strips.
- Prepare the washing buffer by diluting the Wash Buffer 10x (100 mL + 900 mL H<sub>2</sub>O).
- Prepare the antigen by reconstituting the freeze-dried product with the conjugate (volume reported on the label).

Dilute samples 1:101 distributing 10 µL of serum into 1 mL of diluent. Dispense 100 µL of each diluted sample per well (duplicate testing is recommended). Place UNDILUTED calibrators (if possible, in duplicate) in a strip (100 µL in each well). Leave one well for the blank, performed using 100 µL of the substrate mixture.

Wells are covered with protective film and incubated for 45 minutes at 37°C. After washing four times for 30 seconds (300 µl ± 75 µl), add 100 µL of immunocomplex (antigen/monoclonal antibodies labelled with peroxidase) to each well and incubate again for 45 minutes at 37°C, covering the wells with the protective film. The plate is washed again 4 times, as described above. Finally, the substrate is distributed, 100 µL/well.

After 15 minutes at room temperature the enzymatic reaction is stopped with 100 µL of Stop Solution.

The adsorbance (O.D.) is read at 450 nm or 450/620 nm within 30 min.

<b>9. SCHEME OF TEST PROCEDURE</b>
------------------------------------

- |        |   |
|--------|---|
| STEP 1 | Place 100 µL of diluted sample/<br>controls in the wells of the strips. |
|        | -   |
|        | Incubate for 45 min. at 37°C  |
|        | -   |
|        | Wash 4 times (300 µl)   |
|        | -   |
| STEP 2 | Add 100 µL of immunocomplex to each well                                |
|        | -   |
|        | Incubate for 45 min. at 37°C  |
|        | -   |
|        | Wash 4 times (300 µl)   |
|        | -   |
| STEP 3 | Add 100 µL of Substrate to each well                                    |
|        | -   |
|        | Incubate for 15 min. at R.T.  |
|        | -   |
| STEP 4 | Add 100 µL of Stop Solution   |
|        | -   |
|        | Read absorbance at 450 nm within 30 min                                 |

## **10. TEST VALIDATION**

Subtract the value of the blank ( $\leq 0.150$ ) from all the other readings. The O.D. values of the control Cut-off serum when tested in triplicate must be within 25% of the mean value. Disregard any abnormal value and recalculate the mean. The Positive control must have an O.D. at least 1.5 times that of the Cut-Off serum. The ratio between Negative Control and Cut-off must be  $\leq 0.6$ . The O.D. of the cut-off must be  $\geq 0.2$  at 450 nm and  $\geq 0.16$  at 450/620 nm.

## **11. INTERPRETATION OF THE RESULTS**

### **Qualitative results**

If the adsorbance of the sample is higher than that of the Cut-Off, the sample is positive for the presence of specific IgM.

Calculate the ratio between the O.D. value of the sample and that of the Cut-off (INDEX).

The sample is considered:

Positive: if the ratio is  $> 1.2$ .

Doubtful:  $\pm 20\%$  of the Cut-Off.

Negative: if the ratio is  $< 0.8$ .

If the result is doubtful, repeat the test. If it remains doubtful, collect a new serum sample.

## **12. LIMITATIONS OF THE PROCEDURE**

All positive test results require careful interpretation since false positive reactions or heterotypic IgM responses may occur with sera from patients with heterophile-positive mononucleosis, or Varicella Zoster.

A specific IgM response may be observed in reactivation and reinfection as well as with primary infections with CMV. Because of all the complications of serologic diagnosis of congenital infection, virus isolation from urine in the first week of life remains the best way to diagnose intrauterine involvement. Absence of CMV specific IgM does not exclude the possibility of CMV infection. It has been reported that 10-30% of infants may fail to develop CMV IgM antibody response despite congenital infection with CMV.

The test results should be used in conjunction with information available from the clinical evaluation and other diagnostic procedures.

Samples which are strongly positive for the presence of anti-Varicella Zoster (VZV) and Epstein Barr IgM antibodies can give false positive results.

## **13. ANALYTICAL SPECIFICITY**

44 serum samples containing potential interfering substances were tested:

- Rheumatoid factor (n=8)
- Heterophile antibodies (n=4)
- Bilirubin (n=8)
- Triglycerides (n=9)
- Varicella Zoster IgM positive (n=2)
- Hypergammaglobulinemia (n=13).

Two doubtful results were found with the VZV and heterophyl-positive samples, while no interference was noted in the other cases.

## **14. DIAGNOSTIC SENSITIVITY AND SPECIFICITY**

In a clinical trial, 105 samples taken from patients whose clinical picture was compatible with a CMV infection, were analysed with the present kit in parallel with the method taken as reference to establish the sensitivity. With the reference method, 62 samples were effectively positive.

To determine the specificity, 222 samples were used, of which 87 came from donors, 92 were IgG-negative and 43 patients were negative in the reference test.

The results are reported in the following table.

		REFERENCE METHOD	
		+	-
DIESSE	+	61	5
	-	1	217

The Enzywell kit has a sensitivity of 98.4 % and a specificity of 97.8%.

**15. PRECISION****In-run Precision (n=6):**

	MEAN INDEX	ST. DEV.	CV%
Negative control	0.2	0.02	10
CMM2	3.3	0.15	5
CMM3	10.3	0.45	4
Positive 1	3.2	0.27	8
Positive 2	1.6	0.08	5
Positive 3	4.3	0.24	6
Positive 4	2.3	0.19	8
Negative 1	0.2	0.02	10
Negative 2	0.3	0.04	13
Negative 3	0.2	0.02	10
Negative 4	0.4	0.03	8
Negative 5	0.3	0.04	13
Negative 6	0.2	0.03	15
Negative 7	0.3	0.04	13
Negative 8	0.3	0.04	13

**Between-run Precision (n=6):**

	MEAN INDEX	ST. DEV.	CV%
Negative control	0.2	0.01	5
CMM2	3.3	0.01	0
CMM3	10	0.36	4
Positive 1	3.2	0.08	3
Positive 2	1.6	0.11	7
Positive 3	4.2	0.08	2
Positive 4	2.6	0.53	20
Negative 1	0.3	0.06	20
Negative 2	0.2	0.03	15
Negative 3	0.3	0.03	10
Negative 4	0.4	0.07	18
Negative 5	0.3	0.04	13
Negative 6	0.2	0.03	15
Negative 7	0.2	0.03	15
Negative 8	0.3	0.04	13

**16. TROUBLE SHOOTING GUIDE**

PROBLEM	POSSIBLE SOURCE	TEST OR ACTION
Invalid run (all negative)	One or more reagents not added or added in wrong sequence	Recheck procedure Check for unused solutions. Repeat test.
	Unreactive plate	Check the code on the package containing the plate (see package insert paragraph 4 for correct code).
		Check for moisture in unused plate. (Silica gel dessiccant must be pale yellow). Repeat test
Invalid run (all positive)	Contamination of substrate	Take new aliquot of substrate.

	Inadequate washing	Ensure that wash apparatus works well
Poor precision	Incomplete washing of wells	Ensure that wash apparatus works well
	Inadequate aspiration of wells	Ensure that wash apparatus works well
	Pipetting error	Check pipette function
	Reagent addition too slow	Avoid drying of the plate after washing step. Add reagents immediately
	Presence of bubbles	Avoid air bubbles during pipetting.
	Optical pathway not clean	Check instrument light source and detector for dirt. Wipe bottom of plate with soft tissue.
Inadequate Color development	Incorrect incubation times or temperature	Check for temperature control and time monitoring
		Adhere to recommended instruction for use.
	Inadequate volume of substrate added to the plate	Check pipette function.

## **17. REFERENCES**

- 1.G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
- 2.H.O. Kangro: Cytomegalovirus serology: does it give the answer? Serodiagnosis and Immunotherapy 1: 91 (1987).
- 3.Grint P.C.A. et al.: Screening tests for antibodies to cytomegalovirus: an evaluation of five commercial products. J. Clin. Pathol. 38: 1059 (1985).
- 4.Van Loon A.M. et al.: Direct enzyme-linked immunosorbent assay that uses peroxidase-labelled antigen for determination of immunoglobulin M antibody to cytomegalovirus. J. Clin. Microbiol. 13: 416
- 5.M. Musiani et al.: Rapid detection of antibodies against cytomegalovirus induced immediate early and early antigens by an enzyme linked immunosorbent assay. J. Clin. Pathol. 37: 122 (1984).
- 6.F. de Ory et al.: Serological diagnosis of cytomegalovirus infections: comparison of six commercial methods of ELISA. Serodiagnosis and Immunotherapy 2: 423 (1988).
- 7.P. Dal Monte et al.: Cytomegalovirus umano. Diagnosis 2: 67 (1990).
- 8.R. Ziegelmaier et al.: ELISA. la Ricerca Clin. Lab. 10: 83 (1980).

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## **17. BIBLIOGRAFIA**

- 1.G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
- 2.H.O. Kangro: Cytomegalovirus serology: does it give the answer? Serodiagnosis and Immunotherapy 1: 91 (1987).
- 3.Grint P.C.A. et al.: Screening tests for antibodies to cytomegalovirus: an evaluation of five commercial products. J. Clin. Pathol. 38: 1059 (1985).
- 4.Van Loon A.M. et al.: Direct enzyme-linked immunosorbent assay that uses peroxidase-labelled antigen for determination of immunoglobulin M antibody to cytomegalovirus. J. Clin. Microbiol. 13: 416
- 5.M. Musiani et al.: Rapid detection of antibodies against cytomegalovirus induced immediate early and early antigens by an enzyme linked immunosorbent assay. J. Clin. Pathol. 37: 122 (1984).
- 6.F. de Ory et al.: Serological diagnosis of cytomegalovirus infections: comparison of six commercial methods of ELISA. Serodiagnosis and Immunotherapy 2: 423 (1988).
- 7.P. Dal Monte et al.: Cytomegalovirus umano. Diagnosi 2: 67 (1990).
- 8.R. Ziegelmaier et al.: ELISA. la Ricerca Clin. Lab. 10: 83 (1980).

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## **17. REFERÊNCIAS**

1. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. H.O. Kangro: Cytomegalovirus serology: does it give the answer? Serodiagnosis and Immunotherapy 1: 91 (1987).
3. Grint P.C.A. et al.: Screening tests for antibodies to cytomegalovirus: an evaluation of five commercial products. J. Clin. Pathol. 38: 1059 (1985).
4. Van Loon A.M. et al.: Direct enzyme-linked immunosorbent assay that uses peroxidase-labelled antigen for determination of immunoglobulin M antibody to cytomegalovirus. J. Clin. Microbiol. 13: 416
5. M. Musiani et al.: Rapid detection of antibodies against cytomegalovirus induced immediate early and early antigens by an enzyme linked immunosorbent assay. J. Clin. Pathol. 37: 122 (1984).
6. F. de Ory et al.: Serological diagnosis of cytomegalovirus infections: comparison of six commercial methods of ELISA. Serodiagnosis and Immunotherapy 2: 423 (1988).
7. P. Dal Monte et al.: Cytomegalovirus umano. Diagnosis 2: 67 (1990).
8. R. Zieglermaier et al.: ELISA. la Ricerca Clin. Lab. 10: 83 (1980).

	EN ES IT	Date of manufacture Fecha de fabricación Data di fabbricazione	FR GR PT	Date de fabrication Ημερομηνία Παραγωγής Data de fabrico
	EN ES IT	Use By Fecha de caducidad Utilizzare entro	FR GR PT	Utiliser jusque Ημερομηνία λήξης Prazo de validade
	EN ES IT	Caution, consult accompanying documents Atención, ver instrucciones de uso Attenzione, vedere le istruzioni per l'uso	FR GR PT	Attention voir notice d'instructions Προειδοποίηση, συμβουλευτείτε τα συνοδά έντυπα Atenção, consulte a documentação incluída
	EN ES IT	Manufacturer Fabricante Fabbricante	FR GR PT	Fabricant Κατασκευαστής Fabricante
	EN ES IT	Contains sufficient for <n> tests Contenido suficiente para <n> ensayos Contenuto sufficiente per "n" saggi	FR GR PT	Contenu suffisant pour "n" tests Περιεχόμενο επαρκές για «n» εξετάσεις Conteúdo suficiente para "n" ensaios
	EN ES IT	Temperature limitation Límite de temperatura Limiti di temperatura	FR GR PT	Limites de température Περιορισμοί θερμοκρασίας Limites de temperatura
	EN ES IT	Consult Instructions for Use Consulte las instrucciones de uso Consultare le istruzioni per l'uso	FR GR PT	Consulter les instructions d'utilisation Συμβουλευτείτε τις οδηγίες χρήσης Consulte as instruções de utilização
	EN ES IT	Biological risks Riesgo biológico Rischio biologico	FR GR PT	Risques biologiques Βιολογικοί κίνδυνοι Risco biológico
	EN ES IT	Catalogue number Número de catálogo Numero di catalogo	FR GR PT	Référence du catalogue Αριθμός καταλόγου Referência de catálogo
	EN ES IT	In Vitro Diagnostic Medical Device Producto sanitario para diagnóstico in vitro Dispositivo medico-diagnostico in vitro	FR GR PT	Dispositif médical de diagnostic in vitro In Vitro Διαγνωστικό Ιατροτεχνολογικό προϊόν Dispositivo médico para diagnóstico in vitro
	EN ES IT	Batch code Código de lote Codice del lotto	FR GR PT	Code du lot Αριθμός Παρτίδας Código do lote

**CE**  
**0123**

**DIESSE Diagnostica Senese**  
**Via delle Rose 10**  
**53035 Monteriggioni (Siena) Italy**  
**Tel. 0577-587111**