*ENZYWELL*

**HSV 2 SCREEN**

**REF** 91019 (96 tests)

Manufactured by: DIESSE Diagnostica Senese  
Via delle Rose 10  
53035 Monteriggioni (Siena) - Italy

**INDICE / INDEX**

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## **17. BIBLIOGRAFIA**

1. W.E. Rawls et al. Measurement of antibodies to herpesvirus types 1 and 2 in human sera. J. Immunol. 104: 599 (1970).
2. J.J. Gibson et al. A cross-sectional study of herpes simplex virus types 1 and 2 in college students: occurrence and determinants of infection. J. Inf. Dis. 162: 306 (1990).
3. R.M. Coleman et al. Determination of Herpes Simplex Virus antibodies by enzyme-linked immunosorbent assay. J. Clin. Microbiol. 18: 287 (1983).

***DIESSE Diagnostica Senese***  
***Via delle Rose 10***  
***53035 Monteriggioni (Siena) Italy.***  
***Tel. 0577-587111***

**NIESSE**  
  
**INSTRUCTIONS FOR USE**

**ENZYWELL**  
**HSV-2 SCREEN**

**REF** 91019

**1. INTENDED USE:**

**IMMUNOENZYMATIC METHOD FOR THE QUALITATIVE DETERMINATION OF IMMUNOGLOBULINS AGAINST HERPES SIMPLEX VIRUS (TYPE 2.1) IN HUMAN SERUM.**

**2. SUMMARY AND EXPLANATION OF THE TEST:**

The Herpes simplex virus (HSV) is a member of the Herpesviridae family, of which two types are known: type 1 (HSV-1) and type 2 (HSV-2) which present slight antigenic differences. HSV-1 causes chiefly oral-facial lesions, while HSV-2 is mainly responsible for genital lesions, but this distinction is not binding, both types occasionally causing infection in either anatomical site. HSV may also cause a form of ocular cheratitis, and lesions of the central nervous system.

HSV can affect practically the whole population. The primary infection is often in a subclinical form and is rarely diagnosed. After a latency period of variable duration, reactivation may occur and viral replication may or may not give rise to clinical lesions. Infection contracted during birth is of particular interest, this being an important cause of morbidity and mortality. It is therefore important to determine the immunitary state of women during pregnancy in order to detect serum conversion.

**3. PRINCIPLE OF THE TEST:**

The test for the assay of anti-Herpes Simplex Virus (HSV) type 1 immunoglobulins (total Ig) is based on the principle of analysis by competition. The anti-HSV-2 antibodies present in the sample compete with the marker (human IgG anti-HSV-2 conjugated with peroxydase) to occupy the few binding sites available on the antigen coated on the solid phase. The higher the antibody concentration in the sample, the lower the number of conjugated antibodies which are bound. The unbound components are eliminated by washing, and the bound enzymatic activity is determined colorimetrically by transformation of a chromogen substrate (the blue colouring of the reaction turns yellow after the addition of the Stop reagent). The intensity of the colour which develops is inversely proportional to the antibody concentration in the sample.

**4. KIT CONTENTS AND REAGENT PREPARATION:**

- Reagents are sufficient for 96 determinations.
- **Bring reagents to room temperature before use.**

**MT PLATE** MICROPLATE. 6 x 16 wells coated with HSV Type 2.1.

Use: open the package at the opposite end from the code (D followed by the lot number) which is useful for identification purposes, remove the support and wells to be used from the foil package, and place the unused wells in the polythene bag with the silica gel, expel the air and seal by pressing the closure.

**CAL** CALIBRATOR. 1 x 0.8 mL.

Contents: Bovine serum in Phosphate buffer 0.01 mol/L with phenol 0.05% and Bronidox 0.02%, liquid, ready for use.

**CONJ** CONJUGATE. 2 x 8 mL. Ready for use.

Contents: anti-HSV-2 monoclonal antibodies labelled with Peroxydase, in phosphate buffer with phenol 0.05% and Bronidox 0.02%.

**WASH BUF 10x** WASH BUFFER 10X (PF93603). 1 x 100 mL. **INTERCHANGEABLE BETWEEN LOTS**

Contents: Phosphate buffered saline, concentrated 10 times; contains Brij 0.5%.

Preparation: dilute the required volume 1:10 with distilled or deionized water in order to obtain the washing buffer ready for use. If crystals are present, they should be dissolved at 37°C before dilution.

**SUBS TMB** SUBSTRATE (PF93619). 1 x 12 mL. Ready for use. **INTERCHANGEABLE BETWEEN LOTS**  
 Contents: Tetramethylbenzidine 0.26 mg/mL and hydrogen peroxide 0.01% stabilized in citrate buffer 0.05 mol/L (pH 3.8).

**H<sub>2</sub>SO<sub>4</sub> 0.3 M** STOP SOLUTION (PF93602). 1x16 mL. **INTERCHANGEABLE BETWEEN LOTS**.  
 H<sub>2</sub>SO<sub>4</sub> 0.3 mol/L, in solution ready for use.

ADHESIVE FILMS (2).  
 POLYTHENE BAG (1).

#### **MATERIALS REQUIRED BUT NOT PROVIDED.**

- Incubator at 37°C
- Microplate reader, wavelength 450 or 450/620 nm, with OD linearity up to 2,000 (at least).
- Microplate washer (preferable) able to dispense volumes in the range 225-375 µl
- Distilled or deionized water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 10, 100, 1000 µl solution
- Disposable gloves
- Timer
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials
- Absorbent tissue.

#### **5. STORAGE AND STABILITY OF REAGENTS**

Reagents must be stored at 2/8°C.

The expiry date is printed on each component and on the box label.

#### **Reagents have a limited stability after opening and/or preparation**

REAGENT	CONDITIONS
MICROPLATE	6 weeks at 2/8°C in the polythene bag
CALIBRATORE	1 week at room temperature
CONJUGATE	1 week at room temperature.
SUBSTRATE	until the expiry date at 2/8°C, 1 week at 15-30°C in the dark
WASH BUFFER	2 weeks at 2/8°C, 5 days at 15/30°C.
STOP SOLUTION	until the expiry date at 2/8°C

#### **6. PRECAUTIONS AND WARNINGS:**

**FOR IN VITRO DIAGNOSTIC USE. STORE AT 2-8°C.**

#### Health and Safety Information

1. Do not pipette by mouth. Wear disposable gloves and eye protection while handling specimens and performing the assay. Wash hands thoroughly when finished.
2. The following reagents contain low concentrations of harmful or irritant substances:
  - a) The Wash Buffer contains detergents
  - b) The conjugate and controls contain phenol
  - c) The substrate is acid.
 If any of the reagents come into contact with the skin or eyes, wash the area extensively with water.
3. Non-disposable apparatus should be sterilized after use. The preferred method is to autoclave for 1 h at 121°C; disposables should be autoclaved or incinerated.
4. Sulphuric acid required for the Stop Solution and hydrochloric acid 2 M used for washing glassware are corrosive and should be handled with appropriate care. If they come into contact with the skin or eyes, wash thoroughly with water.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1.0%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1.0% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical precautions

1. The extinction value of the calibrator, and also of the samples, may differ slightly from plate to plate. Therefore, if strips taken from different plates, even if from the same batch, are to be used in the same run, the calibrator must be repeated.
2. Allow all reagents and samples to come to room temperature (18-30°C) before use. Immediately after use return reagents to the recommended storage temperature.
3. Do not use the reagents beyond the stated expiry date. Microbiological contamination of reagents must be avoided as this may reduce the life of the product and cause erroneous results.
4. Do not modify the Test Procedure or substitute reagents from other manufacturers or other lots unless the reagent is stipulated as interchangeable. Do not reduce any of the recommended incubation times.
5. Any glassware to be used with the reagents should be thoroughly washed with 2M hydrochloric acid and then rinsed with distilled water or high quality deionized water.
6. Do not expose reagents to strong light or hypochlorite fumes during storage or during incubation steps.
7. Do not allow wells to become dry during the assay procedure.
8. Care must be taken not to cross-contaminate reagents. It is important that pipettes are dedicated for use with the substrate and for use with the conjugate.
9. Care should be taken to avoid touching or splashing the rim of the well with conjugate. Do not "blow-out" from microplates.
10. Enzyme immunoassays can occasionally exhibit an "edge effect" which must be minimized by increasing the humidity during incubation steps. Plates must be covered with their covers and incubated at 37°C either in a water bath with a rack or float to support the plates if necessary, or in an incubator. Alternatively, plates can be incubated in an approved analyzer. See the appropriate operating manual for further details. CO<sub>2</sub> incubators must not be used.
11. Ensure that the bottom of the plate is clean and dry, and that no bubbles are present on the surface of the liquid before reading the plate.
12. Use of highly hemolyzed samples, incompletely clotted sera, plasma samples containing fibrin or samples with microbial contamination may give rise to erroneous results.
13. Avoid contaminating microwells with the dust from disposable gloves.
14. Use of the kit with automatic instruments must be validated by the user.
15. For each instrument used, read the manufacturer's instructions manual carefully to obtain additional information on the following points:
  - installation and particular requisites
  - operating principles, instructions, precautions and risks
  - manufacturer's specifications and instrument performance
  - servicing and maintenance.

**7. TYPE OF SPECIMENS AND STORAGE:**

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice. The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at -20°C, and can be thawed a maximum of 3 times. Defrosted samples must be carefully mixed before performing the test. Heat inactivation can lead to erroneous results. The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

Strongly lipemic, icteric or contaminated samples should be avoided. If a new sample cannot be obtained, such samples should be clarified by filtration (0.45 µm) or centrifugation (3000 rpm x 10').

**The test is not applicable to human plasma.**

**8. TEST PROCEDURE:**PREPARATION

Before testing, bring all reagents and samples to room temperature (18-25°C).

Preparation of Washing buffer: 18 mL of distilled water + 2 mL Wash Buffer 10x per strip.

EXECUTION OF THE TEST

The extinction value of the calibrator, and also of the samples, may differ slightly from plate to plate. Therefore, if strips taken from different plates, even if from the same batch, are to be used in the same test, for a correct evaluation of the samples being tested the calibrator should be repeated.

1. Distribution of the samples:  
Dispense 30 µL of calibrator in a well on the plate (preferably in duplicate), and 30 µL of the samples being tested in the remaining wells.  
Add 100 µL of conjugate to each well. Mix carefully.

2. Incubation:  
Incubate the plate covered with the adhesive foil at 37°C for 90 minutes.
3. Rinsing:  
Remove the adhesive film, aspirate the contents of all the wells and rinse 4 times, filling each well with 300 µL of washing solution. Wait 30" between each wash.
4. Distribution of the substrate:  
Dispense 100 µL of the substrate in each well.
5. Substrate incubation:  
Incubate the plate for 15 minutes at room temperature (18-25°C).
6. Interruption of the reaction:  
Dispense 100 µL of blocking solution, in the same order as that followed for point 4.
7. Reading:  
Take the photometric readings within 30 minutes at 450 nm.

## **9. TEST PROCEDURE FOR HSV-2 SCREEN**

- STEP 1            Place 30 µL of calibrator and serum sample in the wells of the strips and add 100 µL of conjugate.  
Mix well.
- Incubate for 90 min. at 37°C
- Wash 4 times (300 µL)
- STEP 2            Place 100 µL of Substrate in each well
- Incubate for 15 min. at room temperature
- STEP 3            Add 100 µL of Stop Solution
- Read absorbance at 450 nm within 30 min

## **10. VALIDATION OF THE TEST**

O.D. Calibrator  $\geq 0.6$  at 450 nm or  $\geq 0.56$  at 450/620 nm.

## **11. INTERPRETATION OF THE TEST**

The Cut-Off is calculated by multiplying the O.D. of the Calibrator by a factor which is reported on the label.

Example: Calibrator: O.D. = 1.100            **Factor = 0.6**  
Cut-Off = 1.100 x 0.6 = 0.66.

If the absorbance value of the sample is lower than that of the Cut-Off, the sample is positive for the presence of anti-HSV type 1 immunoglobulins.

It is advisable to confirm each positive result by repeating the test. If this second analysis does not confirm the first, the latter is considered valid and the sample is negative. If the second result is doubtful, a new sample should be tested.

## **12. LIMITATIONS**

The test is not able to discriminate between the presence of IgG and that of IgM.

The test cannot be used on its own for clinical diagnostic purposes. A negative result does not exclude the possibility of the disease.

The Test Procedure and Interpretation of the Test must be followed. Non-repeatable reactive results may be obtained with any EIA procedure.

**13. ANALYTICAL SPECIFICITY**

The presence of anti-CMV, anti-EBV, anti-Rubella or anti-Mumps antibodies does not alter the kit's performance.

**14. DIAGNOSTIC SENSITIVITY AND SPECIFICITY**

217 samples were analysed in parallel with other commercial kits. Samples found to be doubtful in all the tests were considered as positive. The following results were obtained:

	OTHER METHODS	
	POS	NEG
DIESSE POS	30	2
NEG	11	174

The diagnostic specificity was 99% (IC at 95% 96-100%) and the diagnostic sensitivity was 73% (IC at 95%, 59-95%).

**15. PRECISION**

Table 1 "In run" Precision

Sample	N	O.D.	CV%
Neg. Calibrator	16	1.146	6
1 (neg)	4	0.778	8
2 (wak pos.)	4	0.249	5
3 (pos.)	4	0.1	4

Tab. 2 Precision between runs and between lots (n=3)

	O.D. Sample/Cut-off			Media	S.D.	CV%
Neg	1.193	1.131	0.143	1.156	0.03	3
Weak pos.	0.485	0.362	0.404	0.417	0.06	15
Pos.	0.205	0.145	0.172	0.174	0.03	17

**16. TROUBLE SHOOTING GUIDE**

PROBLEM	POSSIBLE SOURCE	TEST OR ACTION
Invalid run (all negative)	One or more reagents not added or added in wrong sequence	Recheck procedure Check for unused solutions. Repeat test.
	Unreactive plate	Check the code on the package containing the plate (see package insert point 4 for correct code).
		Check for moisture in unused plate. (Silica gel dessiccant must be pale yellow). Repeat test
Invalid run (all positive)	Contamination of substrate	Take new aliquot of substrate.
	Inadequate washing	Ensure that wash apparatus works well
Poor precision	Incomplete washing of wells	Ensure that wash apparatus works well
	Inadequate aspiration of wells	Ensure that wash apparatus works well
	Pipetting error	Check pipette function
	Reagent addition too slow	Avoid drying of the plate after washing step. Add reagents immediately
	Presence of bubbles	Avoid air bubbles during pipetting.
	Optical pathway not clean	Check instrument light source and detector for dirt. Wipe bottom of plate with soft

		tissue.
Inadequate Color development	Incorrect incubation times or temperature	Check for temperature control and time monitoring
		Adhere to recommended instruction for use.
	Inadequate volume of substrate added to the plate	Check pipette function.

## **17. REFERENCES**

1. W.E. Rawls et al. Measurement of antibodies to herpesvirus types 1 and 2 in human sera. J. Immunol. 104: 599 (1970).
2. J.J. Gibson et al. A cross-sectional study of herpes simplex virus types 1 and 2 in college students: occurrence and determinants of infection. J. Inf. Dis. 162: 306 (1990).
3. R.M. Coleman et al. Determination of Herpes Simplex Virus antibodies by enzyme-linked immunosorbent assay. J. Clin. Microbiol. 18: 287 (1983).

***DIESSE Diagnostica Senese***  
***Via delle Rose 10***  
***53035 Monteriggioni (Siena) Italy.***  
***Tel. 0577-587111***