Instructions for Use

380316-A

English

OC-FCa Reagent

REF V-PH11

INTENDED USE

For quantitative measurement of calprotectin in faeces using automated immunochemical analyser [OC-SENSOR PLEDIA].

INTRODUCTION

Calprotectin is an inflammatory protein abundantly present in neutrophils¹. The amount of calprotectin in faeces increases accompanied with intestinal mucosal inflammation. Therefore measuring the amount of calprotectin in faeces can be used for disease monitoring of IBD patients and distinguishing inflammatory bowel disease (IBD; e.g. ulcerative colitis, Crohn's disease) from functional intestinal disorders (e.g. Irritable bowel syndrome)²⁻⁵. Currently, the methods used for detection of calprotectin include immunoassay latex agglutination, ELISA, FEIA, Chemiluminescent immunoassay and Immunochromatography. OC-FCa Reagent is an immunoassay test reagent used together with an automatic analyser for measurement of calprotectin in faeces. OC-FCa Reagent contains the latex particles coated with anti-human calprotectin antibodies for optical measurement of the latex agglutination reaction.

PRINCIPLE OF THE METHOD

The test method is based on a latex agglutination reaction. A latex reagent is prepared by coating anti-human calprotectin antibodies to polystyrene latex particles. When this reagent is mixed with the sample, the anti-human calprotectin antibodies, which were attached to latex, react with the calprotectin in the sample, and the latex aggregate is formed in the latex agglutination reaction. The change in absorbance per unit time resulting from the latex agglutination reaction is proportional to the concentration of calprotectin in the sample. A dose response curve of the absorbance unit (OD) vs. concentration is generated using the results obtained from the calibrators. The concentration of calprotectin in the patient sample is determined from this curve.

REAGENTS

Reagents are stable until the date printed on the label assuming the container remains unopened at a storage temperature of 2-10°C.

OC-FCa Latex Reagent8mL, 2 vials, maximum 50 tests/vial
(Latex particles suspension coated with anti-human calprotectin antibodies)

OC-FCa Buffer15mL, 2 vial, maximum 50 tests/vial
(N-2-Hydroxyethylpiperazine-N'-2-ethanesulfonic acid (HEPES) buffered saline)

WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use only.
- Before performing measurement, be sure to read the instruction manual of the analyser.
- Fully prepare and adjust the analyser before measurement.
- Create a calibration curve for each measurement. Also be sure to create a new calibration curve when OC-FCa Reagent from a different lot is used.
- Use the Calibrator (REF V-PH12) for the creation of calibration curve.
- If the measurement results exceed the measurement range, use dedicated sample diluent to dilute the sample and perform measurement again.
- Store all reagents at 2-10°C and avoid freezing.
- Do not use reagents that have passed their expiration date.
- The test sample may contain microorganisms. Therefore, use caution when handling. After use, all samples and other materials must be considered as medical waste and disposed of accordingly.
Example of treatment: Soak for 1 hour or longer in a sodium hypochlorite solution (available chlorine concentration 1000 ppm or greater). (Neutralize any substances that contain acids before soaking.) Alternatively, treat in an autoclave at 121°C for 20 minutes. (Do not treat any items which sodium hypochlorite have adhered to in this way.)
- Disposal of used reagents and containers as medical waste in accordance with local regulations.
- If the product is used in any way other than that specified here, the reliability of measurement results cannot be guaranteed. Be sure to follow the procedure.
- A clinical diagnosis based on the measurement results must be a comprehensive judgment made by the attending physician, including factors such as clinical symptoms and other test results.

SAMPLE COLLECTION DEVICE

- Use the following designated sampling container for collecting specimens.
OC-Auto Sampling Bottle 3 (REF V-PZ25, V-PZ26)
- Collect faecal samples by scraping the surface of the stools in different areas. Collect the amount necessary to cover the groove of the probe.
- Check that the test sample has become fully suspended in the buffer solution inside the sampling bottle.
- Received specimens should be stored at 2-8°C and collected within 3 days to designated sampling container⁶.
- Calprotectin in the OC-Auto Sampling Bottle 3 is stable at room temperature for 3 days, at 2-10°C for 14 days⁷.

PREPARATION OF REAGENTS

- OC-FCa Latex Reagent
Ready to use.
Before measurement, bring reagent to room temperature and invert the OC-FCa Latex Reagent gently several times to assure uniform suspension.
After opening, reagent stability (Decreasing DA values $\leq 5\%$):
 - In refrigerator (2-10°C): 14 days
 - Stored onboard : 14 days (Reagent cap needs to be closed when not testing)
 Do not mix reagent with another bottle.
- OC-FCa Buffer
Ready to use.
Before measurement, bring reagent to room temperature and invert the OC-FCa Buffer gently several times.
After opening, buffer stability:
 - In refrigerator (2-10°C): 14 days.
 - Stored onboard : 14 days (Reagent cap needs to be closed when not testing).
 Do not mix buffer with another lot.
Do not mix buffer from different bottles.
Do not use precipitated buffer, precipitate of buffer may cause problems.

REQUIRED MATERIALS, NOT PROVIDED BY THE MANUFACTURER (Prepare these materials before analysis)

- Wash solution: Sodium hypochlorite 0.15% (0.10%-0.30 % is acceptable)
- Purified water (for washing): Distilled or de-ionized water (1.0-10.0M Ω -cm is acceptable)
- Sample cup: To use sample cups, input the specification of the sample cup into the analyser.
- Printer Paper: Use thermal printer paper which fits the analyser.

TEST PROCEDURE

Perform measurements according to the instruction manual of the analyser.

- Set the OC-FCa Latex Reagent and OC-FCa Buffer on board the analyser.
- Input lot number parameters into the analyser.
- Check water, wash solution and drain tank volumes.
- Load control materials and/or sample collection devices.
- Start analysis.
- After measurement is completed, ensure the results of controls are within the acceptable range.

RESULTS

The reaction of each sample is compared to the calibration curve that was created previously. The concentration of calprotectin ($\mu\text{g/g}$ (faeces)) is found and a positive or negative judgment is made.

INTERPRETATION OF RESULTS

Results for the OC-FCa Reagent greater than 50 $\mu\text{g/g}$ (faeces) should be considered positive results. Establish own cut-off values by each laboratory is recommended.

LIMITATIONS

- Calibrator and Controls should be analysed immediately after dispensing into sample cups, in order to avoid evaporation.
- Samples with exceeded 50000 $\mu\text{g/g}$ of calprotectin concentration may have results within measurement range because of prozone effect.
- Taking NSAIDs (non-steroidal anti-inflammatory drug) regularly may raise faecal calprotectin levels⁸.

ANALYTICAL SPECIFICITY

Almost no effect on the measurement value was found from the following co-existing substances at the indicated concentrations.

Substance	Concentration	Substance	Concentration
bilirubin	18mg/dL	Protein(bovine serum albumin)	3.5g/dL
chyle (Formazin turbidity)	1500FTU	Glucose	4.0g/dL
haemoglobin	125mg/dL	barium sulfate	100mg/dL

INTERNAL QUALITY CONTROL

Each laboratory should establish a quality control program to monitor the performance of OC-FCa Latex Reagent. It is recommended to use the following materials for the quality control in your laboratory.

OC-FCa Control LV1 (REF V-PH13)
OC-FCa Control LV2 (REF V-PH14)

PERFORMANCE CHARACTERISTICS

- Sensitivity
When absorbance change of control samples of 0 and 25 $\mu\text{g/g}$ calprotectin was measured, mean - 2.6 standard deviation (SD) of 25 $\mu\text{g/g}$ was greater than mean + 2.6 SD of 0 $\mu\text{g/g}$.
- Accuracy
When control samples of known concentration were measured, the value obtained was within $\pm 15\%$ of the indicated value.
- Inter/Intra-day reproducibility
When the same sample was measured 10 times within the same run, the coefficient of variation (CV) for the values obtained was 5% or less.
When the same control samples (178 $\mu\text{g/g}$ -1373 $\mu\text{g/g}$) were measured for 20days, the coefficient of variation (CV) for the values obtained for inter/intra-day was within 5%.
- Measurement range
20 - 2720 $\mu\text{g/g}$ (faeces)
- Acceptable range between-reagent batches
 - OC-FCa Latex Reagent
When control samples with known concentration were measured, the value obtained was within 100 $\pm 15\%$ of the indicated value.
 - OC-FCa Buffer
When OC-FCa Calibrator was measured, the OD value obtained was within 100 $\pm 5\%$ as compared with the OD value obtained by measuring OC-FCa Calibrator using in-house standard products.

PRODUCT CODE, PRODUCT NAME & STORAGE

Product code	Product name	Contents	Storage
V-PH11	OC-FCa Reagent	2x8mL 2x15mL	2-10°C
V-PH12	OC-FCa Calibrator	6x1mL	2-8°C
V-PH13	OC-FCa Control LV1	2x5mL	2-8°C
V-PH14	OC-FCa Control LV2	2x5mL	2-8°C
V-PZ25	OC-Auto Sampling Bottle 3	100 bottles	1-30°C
V-PZ26	OC-Auto Sampling Bottle 3 without barcode	100 bottles	1-30°C

REFERENCE

- I. Dale et al.: European Journal of Biochemistry, 134: 1-6, 1983.
- J. Tibble et al.: Gut, 47: 506-513, 2000.
- J. Tibble et al.: Gastroenterology, 123: 450-460, 2002.
- M. Montalto et al.: The American Journal of Gastroenterology, 102: 132-136, 2007.
- A. Schoepfer et al.: The American Journal of Gastroenterology, 105: 162-169, 2010.
- H. Tøn et al.: Clinica Chimica Acta, 292, 41-54, 2000.
- S. Hiraoka et al.: Intestinal Research, 17, 202-209, 2018.
- J. Tibble et al.: Gut, 45: 362-366, 1999.

Advena Ltd.

Tower Business Centre, 2nd Flr., Tower Street, Swatar, BKR 4013 Malta

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EIKEN CHEMICAL CO., LTD.